

Patient transfer out of bed

When lateral and seated transfers out of bed are required, a range of patient handling equipment options are available to help the caregiver.

Lateral Transfers

When a transfer out of bed in a supine position is required, e.g. onto a stretcher chair, tilt table or a multi-position aid like Sara Combilizer®, a lateral transfer aid or a patient lift system can be used.



Maxi Slides®

- Sliding sheets made from a strong, ultra low friction material with excellent gliding properties in both directions.
- A choice of single patient use Maxi Slide™ Flites or washable, reusable Maxi Slides™ help caregivers with everyday patient repositioning.



Maxi Sky 2 / Maxi Move®

- Maxi Sky 2 ceiling lift and Maxi Move floor lifter.
- Seated transfers and lateral transfers with a variety of spreader bars and Arjo slings specific to the needs of the patient.
- Maxi Sky 2 can be used for walking and stepping practice with walking slings protecting the patient from falls during such activities.
- Appropriate sling selection matching the assessment of Patient Readiness is fundamental to ensure patients at varying levels are transferred securely.

Maxi Air®

- Single patient use air assisted mattress system for lateral transfers.
- The system includes a perforated inflatable mattress and an air supply pump.
- The perforated underside of the mattress enables air to flow through the mattress to create a cushion that helps lighten static load, reducing the push-pull forces and allowing the patient to be transferred comfortably and securely.



Maxi Transfer Sheet

- Used together with Maxi Sky 2 or Maxi Move.
- The dual purpose Maxi Transfer Sheet is designed to replace the hospital bed sheet.
- Combining the benefits of a transfer sling for lateral transfer, and the functionality of bed linen, with soft breathable fabric construction, Maxi Transfer Sheet can remain in place under the patient in between transfers.

